

Product Datasheet

E. coli recQ protein (orb245697)

Description	Recombinant E.coli ATP-dependent DNA helicase recQ
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	This is His-tag protein
Protein Sequence	AQAEVLNLESGAKQVLQETFGYQQFRPGQEEIIDTVLSGRDCLVVMPTGGGKSL CYQIPALLNLGLTVVVSPLISLMKDQVDQLQANGVAAACLNSTQTREQQLEVM GCRTGQIRLLYIAPERLMLDNFLEHLAHWNPNVLLAVDEAHCSQWGHDFRPEYA ALGQLRQRFPTLPFALTATADDTTTRQDIVRLGLNDPLIQISSFDRPNIRYMLME KFKPLDQLMRYVQEQRGKSGIYCNSRAKVEDTAARLQSKGISAAAYHAGLENN VRADVQEKQFQRDDLQIVVATVAFGMGINKPNVRFVVFHDIPRNIESYYQETGRA GRDGLPAEAMLFYDPADMAWLRRCLEEKPPGQLQDIERHKLNAMGAFAEAQTC RRLLVLLNYFGEGRQEPGNCDCIDPPKQYDGDSTDAQIALSTIGRVNQRFMGY VVEVIRGANNQRIRDYGHDKLVYGMGRDKSHEHWVSVIRQLIHLGLVTQNIQAQ HSALQLTEAARPVLRGESSQLAVPRIVALKPKAMQKSFGGNYDRKLFKLRKLR KSIADENVPYVVFNDATLIEMAEQMPITASEMLSVNGVGMKLERFGKPFMAL IRAHVDGD
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	71.9 kDa
Uniprot ID	P15043
Expiration Date	6 months from date of receipt.