

Product Datasheet

TRPV4 antibody (orb2457)

Description	TRPV4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human TRPV4 (301-400/871 aa)
Target	TRP12/TRPV4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	96 kDa
Uniprot ID	Q9HBA0
Dilution Range	IHC-P=1:100-500, Flow-Cyt=1ug/test, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.