

Product Datasheet

HS3ST1 antibody (Biotin) (orb24580)

Description	HS3ST1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Heparan sulfate glucosamine 3-O-sulfotransferase 1 protein (21-307AA)
Target	HS3ST1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O14792
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.