

Product Datasheet

Human HNRNPA2B1 protein (orb245864)

Description	Recombinant human Heterogeneous nuclear ribonucleoproteins A2/B1
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	This is His-tag protein
Protein Sequence	MEREKEQFRKLFIGGLSFETTEESLRNYYEQWGKLTDCVVMRDPASKRSRGFGF VTFSSMAEVDAAAMAARPHSIDGRVVEPKRAVAREESGKPGAHTVKKLFVGGIK EDTEEHHLRDYFEEYGKIDTIEIITDRQSGKKRGFGFVTFDDHDPVDKIVLQKYHTI NGHNAEVRKALSRQEMQEDLEVAILEVAPVMEEEEEEDMVVEDLDMATRVGATE VVMTTMEEEIMEVEITMILEIITSNLLTTVQ
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	30.4 kDa
Uniprot ID	I6L957
Expiration Date	6 months from date of receipt.