

Product Datasheet

Human RAD23A protein (orb245970)

Description	Recombinant human UV excision repair protein RAD23 homolog A
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	This is His-tag protein
Protein Sequence	MAVTITLKTLLQQTFKIRMEPDETVKVLKEKIEAEKGRDAFPVAGQKLIYAGKILS DDVPIRDYRIDEKNFVVVMVTKTKAGQGTSAPPEASPTAAPESSTSFPPAPTSGM SHPPPAAREDKSPSEESAPTTSPESVSGSVPSSGSSGREEDAASTLVTGSEYETM LTEIMSMGYERERVVAALRASYNPNPHRAVEYLLTGIPGSPEPEHGSVQESQVSEQ PATEAAGENPLEFLRDQPQFQNMQRQVIQQNPALLPALLQQLGQENPQLLQQISR HQQFIQMLNEPPGELADISDVEGEVGAIGEEAPQMNYIQVTPQEKEAIERLKAL GFPESLVIQAYFACEKNENLAANFLLSQNFDDE
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	41.6 kDa
Uniprot ID	P54725
Expiration Date	6 months from date of receipt.