

Product Datasheet

Human VCAM1 protein (orb246031)

Description

Recombinant human Vascular cell adhesion protein 1

Endotoxins

Not test.

Preservatives

If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

Form/Appearance

Liquid or Lyophilized powder

Storage

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Tag

N-terminal 6xHis-tagged

Note

For research use only

Application notes

This is His-tag protein

Protein Sequence

FKIETTPESRYLAQIGDSVSLTCSTTGCESPFFSWRTQIDSPLNGKVTNEGTTSTL
TMNPVSFGNEHSYLCTATCESRKLEKGIQVEIYSFPKDPEIHLSGPLEAGKPITVKC
SVADVYPFDRLEIDLLKGDHLMKSQEFLEDADRKSLETKSLEVTFTPVIEDIGKVL
VCRAKLHIDEMDSVPTVRQAVKELQVYISPKNTVISVNPSTKLQEGGSVTMTCCS
EGLPAPEIFWSKKLDNGNLQHLSGNATLTLIAMRMEDSGIYVCEGVNLIGKNRKE
VELIVQEKPFTEISPGPRIAAQIGDSVMLTCSVMGCESPSFSWRTQIDSPLSGKV
RSEGTNSTLTSPVSFENEHSYLCTVTCGHKKLEKGIQVELYSFPRDPEIEMSGGL
VNGSSVTVSCKVPSVYPLDRLEIELLKGETILENIEFLEDMDKSLENKSLEMTFIP
TIEDTGKALVCQAKLHIDDMEFEPKQRQSTQTLVNVAPRDITVLVSPSSILEEGS
SVNMTCLSQGFPAPKILWSRQLPNGELQPLSENATLTLISTKMEDSGVYLCEGIN
QAGRSRKEVELIIQVTPKDIKLTAFPSSESVKEGDTVIISCTCGNVPETWILKKKAET
GDTVLKSIDGAYTIRKAQLKDAGVYECESKNKVGSQLRSLTLDVQGRENNKDYF
SPE

Purity

Greater than 90% as determined by SDS-PAGE.

Source

Yeast

MW

76.2 kDa

Uniprot ID
P19320
Expiration Date

6 months from date of receipt.