

Product Datasheet

Human NAALADL2 protein (orb246168)

Description

Recombinant human Inactive N-acetylated-alpha-linked acidic

Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	This is His-tag protein
Protein Sequence	GYYVHTNCPSDAPSSGTVDPQLYQEILKTIQAEDIKKSFRNLVQLYKNEDDMEISK KIKTQWTSGLGEDVQFVNYSVLLDLPGPSPSTVTLSSSGQCFFHPNGQPCSEEARK DSSQDLLSYAAYSAGKTLKAEVIDVSYGMADDLKRIRKIKNVTNQIALLKLGKLP LLYKLSSLEKAGFGGVLLYIDPCDLPKTVNPSHDTFMVSLNPGGDPSTPGYPSVD ESFRQRSNLTSLLVQPISAPLVAKLISSPKARTKNEACSSLELPNNEIRVVSMQV QTVTKLKTVTNNVVGFMGLTSPDRYIIVGSHHHTAHSYNGQEWASSTAIITAFIRA LMSKVKRGWRPDRITIVFCSWGGAFTGNIGSYEWGEDFKKVLQKNVVAYISLHSP IRGNSSLYPVASPSLQQLVVEKNNFNCTRRACPETNISSIQIGDADYFINHLGV PIVQFAYEDIKTLEGPSFLSEARFSTRATKIEEMDPSFNLHETITKLSGEVILQIANE PVLPFNALDIALEVQNNLKGDPNTHQLLAMALRLRESAELFQSDMRPANDPK ERAPIRIRMLNDILQDMEKSFLVKQAPPGFYRNILYHLDEKTSRFSILIEAWEHCKP LASNETLQEALSEVLNSINSAQVYFKAGLDVFKSVLDGKN
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	74.6 kDa
Uniprot ID	Q58DX5
Expiration Date	6 months from date of receipt.