

Product Datasheet

Human PGM1 protein (orb246228)

Description	Recombinant human Phosphoglucomutase-1
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Application notes	Full length of His-SUMO-tag and expression region is 1-562aa
Protein Sequence	MVKIVTVKTQAYQDQKPGTSGLRKRVKVFQSSANYAENFIQSIISTVEPAQRQEA TLVVGDDGRFYMKIAIQLIARIAAANGIGRLVIGQNGILSTPAVSCIIRKIKAIIGIIL TASHNPPGGPNGDFGIKFNISNGGPAPEAITDKIFQISKTIIEYAVCPDLKVDLGVLG KQQFDLENKFKPFTVEIVDSVEAYATMLRSIFDFSALKELLSGPNRLKIRIDAMHG VVGPPYVKKILCEELGAPANSVNCVPLEDFGGHHPDPNLTYAADLVETMKSGEH DFGAAFDGDGDRNMILGKHGFFVNPSSVAVIAANIFSIPYFQQTGVRGFAFSM PTSGALDRVASATKIALYETPTGWKFFGNLMDASKLSLCGEESFGTGSDHIREKD GLWAVLAWLSILATRKQSVEDILKDHQYGRNFFTRYDYEEVEAEGANKMMK DLEALMFDRSFVGKQFSANDKVYTVEKADNFEYSDPVDGSIERNQGLRLIFTDG SRIVFRLSGTGSAGATIRLYIDSYEKDVAKINQDPQVMLAPLISIALKVSQQLQERTG RTAPTIVT
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	77.4 kDa
Uniprot ID	P36871
Expiration Date	6 months from date of receipt.