

Product Datasheet

Human DTWD1 protein (orb246332)

Description	Recombinant human DTW domain-containing protein 1
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal GST-tagged
Note	For research use only
Application notes	Partial of the full length of 1-304aa of GST-tag and expression region is 1-109aa
Protein Sequence	MSLNPPIFLKRSEENSSKFVETKQSQTTSIASEDPLQNLCLASQEV LQKAQQSGR SKCLKCGGSRMFYCYTCYVPVENVP IEQIPLVKLPLKIDIIKHPNETDGKSTAI
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	39.2 kDa
Uniprot ID	Q8N5C7
Expiration Date	6 months from date of receipt.