

Product Datasheet

Human PTPRB protein (orb246420)

Description	Recombinant human Receptor-type tyrosine-protein phosphatase beta
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Application notes	Full length of Cytoplasmic of His-SUMO-tag and expression region is 1643-1997aa
Protein Sequence	RQKVSHGRERPSARLSIRDRPLSVHLNLGQKGNRKTSCPIKINQFEGHFMKLQA DSNYLLSKEYEELKDVGRNQSCDIALLPENRGKNRYNNILPYDATRVKLSNVDDD PCSDYINASYIPGNNFRREYIVTQGPLPGTKDDFWKMWVEQNVHNIVMTQCVE KGRVKCDHYWPADQDSLYYGDLILQMLSESVLPEWTIREFKICGEEQLDAHRLIR HFHYTVWPDHGVPETTQSLIQFVRTVRDYINRSPGAGPTVVHCSAGVGRTGTFI ALDRILQQLDSKDSVDIYGAVHDLRLHRVHVMVQTECQYVYLHQCVRDVLARKL RSEQENPLFIYENVNPEYHRDPVYSRH
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	57.4 kDa
Uniprot ID	P23467
Expiration Date	6 months from date of receipt.