

Product Datasheet

Human RBMY1F protein (orb246552)

Description	Recombinant human RNA-binding motif protein, Y chromosome, family 1 member F/J
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Application notes	This is His-SUMO-tag protein
Protein Sequence	MVEADHPGKLFIGGLNRETNEKMLKAVFGKHGPISEVLLIKDRTSKSRGFATFEE NPADAKNAAKDMNGTSLHGKAIKVEQAKKPSFQSGGRRRPPASSRNRSPSGSL RSARGSSGGTRGWLPSEHGLDDGGYTPDLKMSYSRGLIPVKRGPSSRSGGPP PKKSAPSAVARNSNSWMGSQGPMQSRRENYGVPPRRATISSWRNDRMSTRHDG YATNDGNHPSCQETRDYAPPSRGYAYRDNGHSNRDEHSSRGYRNHRSSRETR DYAPPSRGHAYRDYGHSSRRDESYSGYRNHRSSRETRDYAPPSRGHGYRDYGH SRRHESYSGYRNHPSSRETRDYAPPHRDYAYRDYGHSSWDEHSSRGYSYHD GYGEALGRDHSEHLGSSYRDALQRYGTSHGAPPARGPRMSYGGSTCHAYSNT RDRYGRSWESYSSCGDFHYCDREHVCRKDQRNPPSLGRVLPDPREAYGSSSYV ASIVDGGESRSEKGDSSRY
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	71.7 kDa
Uniprot ID	Q15415
Expiration Date	6 months from date of receipt.