

Product Datasheet

E.coli poly protein (orb246712)

Description	Recombinant Hepatitis C virus polyprotein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	This is His-tag protein
Protein Sequence	APITAYAQQTRGLFSTIVTSLTGRDNTNENRGEVQVLSTATQSFLGTAVNGVMWT VYHGAGSKTISGPKGPVNQMYTNVDQDLVGWPAPPGVRSLAPCTCGASDLYLV TRHADVVPVRRRGDTRGALLSPRIATLKGS GG PLLCPMGHAAGLFRAAVCTR GVAKAVDFVSVESLETTMRSPVFTDNSTPPAVPQTYQVAHLHAPTGS GKSTKVP AAYAAQGYKVLVLNPSVAATLGFGAYMSKAYGIDPNIRTGVRTITTGAPITYSTYG KFLADGGCGSGGAYDIIICDECHSTDSTILGIGTVLDQAETAGVRLTILATATPPGS VTTPHSNIEEVALPTTGEIPFYGKAIPLELIKGRHLIFCHSKKKCDELAKQLTSLGL NAVAYYRGLDVSVIPTSGDVVVCATDALMTGYTGDFDSVIDCNTSVIQTVD FSL DPTFSIETTTVPQDAVSRSQRRGRTGRGRLGIYRYVTPGERPSGIFDTSVLCECY DAGCAWYELTPAETTTRLRAYFNTPGLPVCQDHLEFWEVFTGLTQIDGHFLSQ TKQTGENFPYLVAYQATVCARALAPPSWDTMWKCLIRLKPTLHGPTPLLYRLGS VQNEVTLTHPITKYIMACMSADLEVVT
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	71.3 kDa
Uniprot ID	S4UAW6
Expiration Date	6 months from date of receipt.