

Product Datasheet

Human ARNT protein (orb246723)

Description	Recombinant Homo sapiens Aryl hydrocarbon receptor nuclear translocator
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	This is His-tag protein
Protein Sequence	MAATTANPEMTSDVPSLGPASGNSGPGIQGGGAIVQRAIKRRPGLDFDDDDGE GNSKFLRCDDDDQMSNDKERFARSDDEQSSADKERLARENHSEIERRRRNKMTA YITELSDMVPTCSALARKPDKLTILMAVSHMKSLRGNTSTDGSKPSFLTDQ ELKHLILEAADGFLFIVSCETGRVVYVSDSVTPVLNQPQSEWFGSTLYDQVHPDD VDKLREQLSTSENALTGRILDKTGTVKKEGQQSSMRMCMGSRFSICRMRCGS SSVDPVSVNRLSFVRNRCRNLGSKVDGEPHFVVHCTGYIKAWPPAGVSLPD DDPEAGQGSKFCLVAIGRLQVTSSPNCTDMSNVCQPTFISRHNIEGIFTFVDHR CVATVGYQPQELLGKNIVEFCHPEDQQLLRDSFQQVVKLGQVLSVMFRFRSKN QEWLWMRTSSFTFQNPYSDEIEYIICTNTNVKNSSQEPRPT
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	54.8 kDa
Uniprot ID	P27540
Expiration Date	6 months from date of receipt.