

Product Datasheet

Human MYH6 protein (orb246766)

Description	Recombinant Homo sapiens Myosin-6
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	This is His-tag protein
Protein Sequence	NAYQYMLTDRENQSILITGESGAGKTVNTRKVIQYFASIAAIGDRGKKDNANANK GTLEDQIIQANPALEAFGNAKTVRNDNSSRFGKFIRIHFGATGKLASADIETYLLEK SRVIFQLKAERNYHIFYQILSNKKPELLDMLLVNPNPYDYAFVSQGEVSVASIDDS EELMATDSAFDVLGFTSEEKAGVYKLTGAIMHYGNMKFKQKQREEQAEPDGTE DADKSAYLMGLNSADLLKGLCHPRVKVGNEYVTKGQSVQQVYYSIGALAKAVYE KMFNWMVTRINATLETQPRQYFIGVLDIAGFEIFDFNSFEQLCINFTEKLLQOFF NHHMFVLEQEEYKKEGIEWTFIDFGMDLQACIDLIEKPMGIMSILEEECMFPKATD MTFKAKLYDNHLGKSNNFQKPRNIKKGQEAHFSLIHYAGTVDYNILGWLEKNKD PLNETVVALYQKSSLKMATLFSSYATADTGDSGSKGGKKGSSFTQVSALHR ENLNKLMTNLRTTHPHFVRCIIPNERKAPGVMDNPLVMHQLRCNGVLEGIRICRK GFPNRILYGDFRQRYRILNPVAIPEGQFIDSRKGTEKLLSSLDIDHNQYKFGHTKV FFKAGLLGLEEMRDERLSRIITRMQAQARGQLMRIEFKKIVERRDALLVI
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	76.8 kDa
Uniprot ID	P13533
Expiration Date	6 months from date of receipt.