

Product Datasheet

Irf5 antibody (orb24678)

Description	Irf5 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Mouse Interferon regulatory factor 5 protein (1-497AA)
Target	Irf5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P56477
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.