

Product Datasheet

Mouse Mlkl protein (orb246878)

Description	Recombinant Mus musculus Mixed lineage kinase domain-like protein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	This is His-tag protein
Protein Sequence	MDKLGQIIKLGQLIYEQCEKMKYCRKQCQRLGNRVHGLLQPLQRLQAQGKKNLP DDITAAALGRFDEVLKEANQQIEKFSKKSHIWKFSVSGNDKILFHEVNEKLRDVWE ELLLLLQVYHWNTVSDVSQPASWQQEDRQDAEEDGNENMKVILMLQLQISVEEI NKTLKQCSLKPTQEIPQDLQIKEIPKEHLGPPWTKLKTSKMSTIYRGEYHRSPVTIK VFNNPQAESVGIVRFTFNDEIKTMKKFDSNLRIFGICIDQTVKPPEFSIVMEYCEL GTLRELLDREKDLTMSVRSLLVLRAARGLYRLHHSETLHRNIISSSSFLVAGGYQV KLAGFELSKTQNSISRTAKSTKAERSSSTIYVSPERLKNPFCLYDIKAEIYSFGIVLW EIATGKIPFEGCDSSKIRELVAEDKKQEPVGQDCPELLREIINECRAHEPSQRPSV DGRSLSGRERILERLSAVEESTDKKV
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	56.3 kDa
Uniprot ID	Q9D2Y4
Expiration Date	6 months from date of receipt.