

Product Datasheet

Human TUFM protein (orb246889)

Description	Recombinant Homo sapiens Elongation factor Tu, mitochondrial
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Application notes	Full length of His-SUMO-tag and expression region is 44-452aa
Protein Sequence	AVEAKKTYVRDKPHVNVGTIGHVDHGGKTTLTAAITKILAEGGGAKFKKYEEIDNAP EERARGITINAAHVEYSTAARHYAHTDCPGHADYVKNMITGTAPLDGCILVVAAN DGPMPQTREHLLARQIGVEHVVVYVKNKADAVQDSEMVELVEIRELLTEFGYK GEETPVIVGSALCALEGRDPELGLKSVQKLLDAVDTYIPVPARDLEKPFLLPVEAV YSVPGRGTVVTGTLERGILKKGDECELLGHSKNIRTVVTGIEMFHKSLERAEAGD NLGALVRGLKREDLRRGLVMVKPGSIKPHQKVEAQVYILSKEEGGRHKPFVSHF MPVMFSLTWDMACRIILPPEKELAMPGEDLKFNILRQPMILEKGQRFTLRDGNR TIGTGLVTNTLAMTEEEKNIKWG
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	61 kDa
Uniprot ID	P49411
Expiration Date	6 months from date of receipt.