

Product Datasheet

Mouse Gp2 protein (orb246900)

Description	Recombinant mouse Pancreatic secretory granule membrane major glycoprotein GP2
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Application notes	This is His-SUMO-tag protein
Protein Sequence	STIHQGYGRPRNSSNLDLDCGSPDSPSSGICFDPCQNHTVLNDPTRSTENNDSS VAWCDDNLHGWYRFVGDGGVKMPETCVSVFRCHTSAPMWLSGSHPI LGD GIV SHTACANWNENCCFWRSEVQVKACSEELGEYHVYKLQGTPECSLRYCTDPSTA PKNCEITCRPEEECVFQNNNWSCVCRQDLHVSDSQSLQPLLD CGDNEIKVKLDK CLLGGMGFKEEIIAYLNDRNCNGTMQDEPNWVSM TSPVVANYCGNILEKNGT HAIYRNTLSLATDFIIRDFRVNVNFQ CAYPLDMSVSLETALQPIVSSLTVDVDGAG EFNVKMALFQDQSYTNPYEGA EVLLPVESILYVGVLNLRGDTSRFKLLLTN CYATP SEDRHDPVKYFIKNRCPNQ RDSTINVRENGVSSSRFSVQMFMFAGNYDLVFL HCEVYLC DSTTEQCQPS CSTNRLRSSRPAIDYNRVLDLGPITKRS AQSSATSKGT PHTTGFL LAWPMFFLPVFLALLF
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	73 kDa
Uniprot ID	Q9D733
Expiration Date	6 months from date of receipt.