

Product Datasheet

Human DHX9 protein (orb246920)

Description	Recombinant Human ATP-dependent RNA helicase A
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	C-terminal Flag-Myc-tagged
Note	For research use only
Application notes	Partial of HIS/DDK and expression region is 325-840aa
Protein Sequence	NQVGVPWSPQSNWNPWTSSNIDEGPLAFATPEQISMDLKNELMYQLEQDH DLQAILQERELLVKKFESEILEAISQNSVVIIRGATGCGKTTQVPQFILDFFIQND RAAECNIVVTQPRRISAVSVAERVAFERGEEPGKSCGYSVRFESILPRPHASIMFC TVGVLLRKLEAGIRGISHVIVDEIHERDINTDFLLVVLRDVVQAYPEVRIVLMSATI DTSMFCEYFFNCPIIEVYGRTPVQYFLEDICQMTHFVPPPKDKKKKDKDDGG EDDDANCNLICGDEYGPETRLSMSQLNEKETPFELIEALLKYIETLNVPGAVLVFL PGWNLITYTMQKHLEMNPHFGSHRYQILPLHSQIPREEQRKVFDPPVPGVTKVILS TNIAETSITINDVVYVIDSCKQKVKLFTAHHNMNTNYATVWASKTNLEQRKGRAGR VRPGFCFHLCSRARFERLETHMTPEMFRTPLHEIALSIKLLRLGGIGQFLAKAIEPP PLDAVIEAEHTLRELD
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Mammalian cell
MW	61.6 kDa
Uniprot ID	Q08211
Expiration Date	6 months from date of receipt.