



Product Datasheet

SRRT antibody (orb247035)


Description	Goat polyclonal antibody to SRRT
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	SRRT / ARS2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize both reported isoformsWB: Approx 140kDa band observed in lysates of cell line Daudi and in lysates of Mouse Spleen (calculated MW of 101kDa according to NP_056992.4). The observed molecular weight is explained by the acetylation of the protein and it corresponds to earlier finding
Protein Sequence	DPRAIVEYRDLDPD
Clonality	Polyclonal
MW	101; 101; 100; 100
Entrez	51593
Dilution Range	ELISA: 1:128000, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of staining of Daudi cell lysates



Western blot analysis of staining of Daudi cell lysates



Immunohistochemical staining of Human Uterine Tissue (HUT)