



Product Datasheet

GP1BA antibody (orb247044)


Description	GP1BA antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human GP1BA. Please contact us for the exact immunogen sequence. The peptide is available as orb374852.
Target	GP1BA
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for short term. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified
MW	70 kDa
NCBI	21 , 010326 , 034456
Entrez	2811
Dilution Range	IHC-P: 1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat sple...



IHC-P staining of human breast cancer tissue...



Immunohistochemical staining of pig liver...