

Product Datasheet

OPG antibody (orb247239)

Description

OPG antibody

Species/Host

Rabbit

Reactivity

Bovine, Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, ICC, IF, IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human OPG. Please contact us for the exact immunogen sequence. The peptide is available as orb12899.

Target

OPG

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for short term. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by Protein A affinity chromatography.

MW

43 kDa

Uniprot ID

[O00300](#)

NCBI

[002546](#), [002537](#), [31](#)

Entrez


[4982](#)

Dilution Range


WB: 1:100-500, IHC-P: 1:100-500 (based on 0.5 mg/ml)

Expiration Date


12 months from date of receipt.



Western blot analysis of Rat heart (Lane...)



Immunohistochemical staining of paraffin...



IF image of mouse ovary tissue using OPG...