


Product Datasheet

MDR1 antibody (FITC) (orb247246)

Description	MDR1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	ICC, IF
Target	MDR1 (FITC)
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	141 kDa
Uniprot ID	P08183 , P06795 , P43245
NCBI	NM_000927.4 , NP_000918.2 , NP_000918.21 , NM_000927.41
Entrez	5243
Dilution Range	IF/ICC: 1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of mouse kidney tissue.