


Product Datasheet

CXCL11 antibody (orb247297)

Description	Rabbit polyclonal antibody to CXCL11
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human CXCL11(NP_005400.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	10kDa
Uniprot ID	O14625
Dilution Range	WB: 1:500-2000, IHC-P: 1:100-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human pr...