



Product Datasheet

BLZF1 antibody (orb247433)


Description	Rabbit polyclonal antibody to BLZF1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human BLZF1 (NP_003657.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	45kDa
Uniprot ID	Q9H2G9
Entrez	8548
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of various cell lines and mouse skeletal muscle.



Immunohistochemical staining of human gastric tissue.



Immunofluorescence analysis of U2OS cells.