



Product Datasheet

GNL3 antibody (orb247452)


Description	Rabbit polyclonal antibody to GNL3
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IP, WB
Immunogen	Recombinant fusion protein of human GNL3(NP_996562.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	62kDa
Uniprot ID	Q9BVP2
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:10-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of 293T cells.



Immunofluorescence analysis of A549 cells.



Immunoprecipitation analysis of 200ug of cell extract.