



Product Datasheet

INTS6 antibody (orb247538)


Description	Rabbit polyclonal antibody to INTS6
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human INTS6(NP_036273.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	117KDa
Uniprot ID	Q9UL03
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human mammary gland tissue.



Western blot analysis of extracts of various cell lines.



Immunohistochemical staining of mouse lung tissue.