



Product Datasheet

NFKB1 antibody (orb247557)


Description	Rabbit polyclonal antibody to NFKB1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human NFKB1(NP_001158884.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50KDa/120KDa
Uniprot ID	P19838
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse lu...



Immunofluorescence analysis of U2OS cell...



Western blot analysis of extracts of var...