



Product Datasheet

TALDO1 antibody (orb247614)

Description	Rabbit polyclonal antibody to TALDO1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human TALDO1 (NP_006746.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	38kDa
Uniprot ID	P37837
Entrez	6888
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:10-100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cells...



Western blot analysis of extracts of various tissues...