



Product Datasheet

RBBP5 antibody (orb247713)


Description	Rabbit polyclonal antibody to RBBP5
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein of human RBBP5 (NP_005048.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	70kDa
Uniprot ID	Q15291
Entrez	5929
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of various cell lines and tissues.



Immunohistochemical staining of rat intestine tissue.



Immunohistochemical staining of human colon tissue.