


Product Datasheet


NAA10 antibody (orb247736)

Description	Rabbit polyclonal antibody to NAA10
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human NAA10(NP_003482.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	26kDa
Uniprot ID	P41227
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-100
Expiration Date	12 months from date of receipt.




WB analysis of extracts of various cell

...



Immunofluorescence analysis of U2OS cell...



Immunohistochemistry of paraffin-embedde...