



Product Datasheet

NMT2 antibody (orb247761)


Description	Rabbit polyclonal antibody to NMT2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human NMT2 (NP_004799.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	57kDa
Uniprot ID	O60551
Entrez	9397
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-100
Expiration Date	12 months from date of receipt.



WB analysis of extracts of various cell ...



Immunohistochemical staining of mouse ki...



Immunofluorescence analysis of U2OS cell...