



Product Datasheet

QKI antibody (orb247762)


Description	Rabbit polyclonal antibody to QKI
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human QKI(NP_006766.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	38KDa
Uniprot ID	Q96PU8
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



WB analysis of extracts of various cell ...



Western blot analysis of extracts of var...



Immunofluorescence analysis of MCF7 cell...