

Product Datasheet

KRT17 antibody (HRP) (orb24778)

Description	KRT17 antibody (HRP)
Species/Host	Rabbit
Reactivity	Bovine
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Bovine Keratin, type I cytoskeletal 17 protein (1-441AA)
Target	KRT17
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A1L595
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.