



Product Datasheet

ACTL6B antibody (orb247810)


Description	Rabbit polyclonal antibody to ACTL6B
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein of human ACTL6B (NP_057272.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	47-55kDa
Uniprot ID	O94805
Entrez	51412
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of various cell lines and tissues.



Immunofluorescence analysis of A549 cells.



Immunohistochemical staining of human kidney tissue.