



## Product Datasheet

### POLD3 antibody (orb247870)


|                            |                                                         |
|----------------------------|---------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to POLD3                     |
| <b>Reactivity</b>          | Human, Mouse                                            |
| <b>Conjugation</b>         | Unconjugated                                            |
| <b>Tested Applications</b> | ICC, IF, IHC, WB                                        |
| <b>Immunogen</b>           | Recombinant fusion protein of human POLD3(NP_006582.1). |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.      |
| <b>Concentration</b>       | 1 mg/ml                                                 |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.             |
| <b>Note</b>                | For research use only                                   |
| <b>Isotype</b>             | IgG                                                     |
| <b>Clonality</b>           | Polyclonal                                              |
| <b>MW</b>                  | 51kDa                                                   |
| <b>Uniprot ID</b>          | <a href="#">Q15054</a>                                  |
| <b>Dilution Range</b>      | WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-100       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                         |



Immunofluorescence analysis of U2OS cell...



Western blot analysis of extracts of var...



Immunohistochemical staining of human oo...