



Product Datasheet

AKT antibody (orb247874)


Description	Rabbit polyclonal antibody to AKT
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic peptide of human AKT1(NP_005154.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	57kDa
Uniprot ID	P31749
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human gastric tissue.



Western blot analysis of extracts of various tissues.



Immunohistochemistry of paraffin-embedded tissue sections.