



Product Datasheet

XRCC6 antibody (orb247902)


Description	Rabbit polyclonal antibody to XRCC6
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	A synthetic peptide of human Ku70(NP_001460.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70kDa
Uniprot ID	P12956
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



WB analysis of extracts of various cell ...



Immunofluorescence analysis of A549 cell...



Immunoprecipitation analysis of 100ug ex...