


Product Datasheet

PDGFRA (phospho-Y754) antibody (orb248094)

Description	Rabbit polyclonal antibody to PDGFRA
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y754 of human PDGFRA(NP_006197.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	180kDa
Uniprot ID	P16234
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of Hep...