

Product Datasheet

L2HGDH antibody (orb24810)

Description	L2HGDH antibody
Species/Host	Rabbit
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Bovine L-2-hydroxyglutarate dehydrogenase, mitochondrial protein (53-463AA)
Target	L2HGDH
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A7MBI3
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.