

Product Datasheet

Rabbit IgG Isotype Control Polyclonal Antibody APC Conjugate (orb248102)

Description Isotype control antibodies are negative controls used to measure non-specific binding or background ...

Species/Host Rabbit

Conjugation APC

Tested Applications FACS, IF, IHC

Immunogen This is purified naive rabbit sera and as such there was no immunogen.

Preservatives 1X PBS, pH 7.4

Storage Store the Rabbit IgG isotype control antibody at 4°C. Do not freeze; protect from light.

Note For research use only

Application notes Generally to be used at the same concentration as the sample being controlled for.

Formula 0.5 mg/ml antibody in PBS with 0.05% sodium azide

Isotype Rabbit IgG


Clonality Polyclonal

Purity Protein A Affinity Chromatography


Hazard Information This Rabbit IgG isotype control antibody is available for research use only.

Dilution Range FACS: assay dependent, IHC: assay dependent, IF: assay dependent

Expiration Date 12 months from date of receipt.



FACS staining of human samples using Rab...



FACS staining of mouse samples using Rab...