

Product Datasheet

Beta Tubulin Antibody Loading Control (orb248111)

Description

There are five mammalian members of the tubulin

Hete...

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Primate, Rat

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

A recombinant protein fragment from the C-terminal region of human beta tubulin was used as the immunogen for this antibody.

Preservatives

1X PBS, pH 7.4

Storage

Aliquot the loading control antibody and store frozen at -20°C or colder to avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Provided assay concentrations are suggestions only, antibody titration may be required for optimal results.

Formula

0.5 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Rabbit IgG

Clonality

Polyclonal

Purity

Protein A Affinity Chromatography

Hazard

Information

This beta Tubulin antibody loading control is available for research use only.

Entrez


[203068](#)

Dilution Range


Western blot: 2-4ug/ml

Expiration Date

12 months from date of receipt.



Western blot analysis of human samples u...



Western blot analysis of mouse samples u...