



Product Datasheet

AGR2 Antibody (orb248116)

Description	AGR2 (Anterior gradient protein 2 homolog) is thought to promote breast tumorigenesis and is overexp...
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	A recombinant protein fragment from the N-terminal region of human AGR2 was used as the immunogen for this antibody.
Preservatives	1X PBS, pH 7.4
Storage	Aliquot the AGR2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Provided assay concentrations are suggestions only, AGR2 antibody titration may be required for optimal results.
Formula	0.5 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Hazard Information	This AGR2 antibody is available for research use only.
Entrez	10551
Dilution Range	Western blot: 4-8ug/ml,IHC (FFPE): 3-5ug/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of human samples
U...