


Product Datasheet

CLEC4E Antibody (orb248117)

Description	CLEC4 is a member of the C-type lectin family of proteins that plays a role in inflammatory and immu...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A recombinant protein fragment from the C-terminal portion of human CLEC4E was used as the immunogen for this antibody.
Preservatives	1X PBS, pH 7.4
Storage	Aliquot the CLEC4E antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Provided assay concentrations are suggestions only, CLEC4E antibody titration may be required for optimal results.
Formula	0.5 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG
Clonality	Polyclonal
Purity	Protein A affinity chromatography
Hazard Information	This CLEC4E antibody is available for research use only.
Entrez	26253
Dilution Range	Western blot: 2-5ug/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human samples u...