


Product Datasheet

TRH Receptor Antibody (orb248119)

Description	Thyrotropin-releasing hormone (TRH) is produced by tissues of the nervous system and to a lesser deg...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A recombinant protein fragment from the C-terminal region of TRH receptor was used as the immunogen for this antibody.
Preservatives	1X PBS, pH 7.4
Storage	Aliquot the TRH receptor antibody and store frozen at -20°C or colder to avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Provided assay concentrations are suggestions only, TRH receptor antibody titration may be required for optimal results.
Formula	0.5 mg/ml antibody in PBS with 0.1 mg/ml BSA (US source) and 0.05% sodium azide
Isotype	Rabbit IgG
Clonality	Polyclonal
Purity	Protein A purified antibody
Hazard Information	This TRH Receptor antibody is available for research use only.
Entrez	7201
Dilution Range	Western blot: 1-3ug/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human samples u...