

## Product Datasheet

### Acid Phosphatase Antibody (Isoform 2 specific) (orb248124)

## Description

Acid phosphatases are a broad family of enzymes

ho...

## Species/Host

Rabbit

## Reactivity

Human

## Conjugation

Unconjugated

## Tested

WB

## Applications

## Immunogen

The amino acid sequence used as immunogen for this antibody is specific to isoforms 2 (or B) and E of Acid Phosphatase.

## Preservatives

1X PBS, pH 7.4

## Storage

Aliquot the Acid Phosphatase antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

## Note

For research use only

## Application notes

Provided assay concentrations are suggestions only, Acid Phosphatase antibody titration may be required for optimal results. 1. This antibody is specific to isoforms 2 (or B) and E of Acid Phosphatase.

## Formula

0.5 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

## Isotype

Rabbit IgG

## Clonality

Polyclonal

## Purity

Protein A affinity chromatography

## Hazard

## Information

This Acid Phosphatase antibody is available for research use only.

## Entrez


[52](#)

## Dilution Range

Western blot: 1-3ug/ml

## Expiration Date

12 months from date of receipt.



Western blot analysis of human samples u...