

Product Datasheet

Acid Phosphatase / ACP1 Antibody (orb248125)

Description

Acid phosphatases (AP) are a broad family of

homol...

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

A recombinant protein fragment from the N-terminal portion of Acid Phosphatase / ACP1 was used as the immunogen for this antibody. This AA sequence is common to the A, C and D isoforms.

Preservatives

1X PBS, pH 7.4

Storage

Aliquot the ACP1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Provided assay concentrations are suggestions only, ACP1 antibody titration may be required for optimal results. 1. This antibody will detect the A, C and D isoforms of Acid Phosphatase.

Formula

0.5 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Rabbit IgG

Clonality

Polyclonal

Purity

Protein A affinity chromatography

Hazard

Information

This Acid Phosphatase / ACP1 antibody is available for research use only.

Entrez


[52](#)

Dilution Range

Western blot: 1-3ug/ml

Expiration Date

12 months from date of receipt.



Western blot analysis of human samples u...