


Product Datasheet

DC-SIGN Antibody (orb248128)

Description	DC-SIGN (Dendritic Cell-Specific Intercellular adhesion molecule-3-Grabbing Non-integrin) also known...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A recombinant protein fragment from human DC-SIGN was used as the immunogen for this antibody.
Preservatives	1X PBS, pH 7.4
Storage	Aliquot the DC-SIGN antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Provided assay concentrations are suggestions only, DC-SIGN antibody titration may be required for optimal results.
Formula	0.5 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG
Clonality	Polyclonal
Purity	Protein A purified antibody
Hazard Information	This DC-SIGN antibody is available for research use only.
Entrez	30835
Dilution Range	Western blot: 2-5ug/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human samples u...