


## Product Datasheet

### Viperin Antibody (orb248143)

|                            |                                                                                                                         |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Viral infection induces the interferon (IFN) response which triggers the production of an array of a...                 |
| <b>Species/Host</b>        | Rabbit                                                                                                                  |
| <b>Reactivity</b>          | Canine, Human, Mouse, Rat                                                                                               |
| <b>Conjugation</b>         | Unconjugated                                                                                                            |
| <b>Tested Applications</b> | WB                                                                                                                      |
| <b>Immunogen</b>           | A recombinant protein fragment from the N-terminal region of human Viperin was used as the immunogen for this antibody. |
| <b>Preservatives</b>       | 1X PBS, pH 7.4                                                                                                          |
| <b>Storage</b>             | Aliquot the Viperin antibody and store frozen at -20°C or colder to avoid repeated freeze-thaw cycle.                   |
| <b>Note</b>                | For research use only                                                                                                   |
| <b>Application notes</b>   | Provided assay concentrations are suggestions only, Viperin antibody titration may be required for optimal results.     |
| <b>Formula</b>             | 0.5 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide                                                           |
| <b>Isotype</b>             | Rabbit IgG                                                                                                              |
| <b>Clonality</b>           | Polyclonal                                                                                                              |
| <b>Purity</b>              | Protein A affinity chromatography purified antibody                                                                     |
| <b>Hazard Information</b>  | This Viperin antibody is available for research use only.                                                               |
| <b>Entrez</b>              | <b>91543</b>                                                                                                            |
| <b>Dilution Range</b>      | Western blot: 3-6ug/ml                                                                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                         |



Western blot analysis of human samples u...