


Product Datasheet

CRHR2 Antibody (orb248164)

Description

CHCR2, or Corticotropin releasing hormone receptor 2, is a

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A recombinant protein fragment from the N-terminal portion of human CRHR2 was used as the immunogen for this antibody (1).
Preservatives	1X PBS, pH 7.4
Storage	Aliquot the CRHR2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Provided assay concentrations are suggestions only, CRHR2 antibody titration may be required for optimal results.1. This antibody will detect the alpha (411 AA), beta (438 AA) and gamma (397 AA) isoforms of CRHR2.
Formula	0.5 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG
Clonality	Polyclonal
Purity	Protein A affinity chromatography
Uniprot ID	Q13324
Hazard Information	This CRHR2 antibody is available for research use only.
Entrez	1395
Dilution Range	Western blot: 0.5-2ug/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human samples u...